



Product Catalogue.



EUROSILICONE®

NAGÔR®

Providing choices for every patient

GC Aesthetics™ is a leading medical aesthetics company with over three decades of experience manufacturing silicone implants through our brands Nagor® and Eurosilicone®. Thousands of patients worldwide have used our products.

We work in partnership with the medical community to continually develop high quality products and services that meet both the surgeon's and the patient's needs. As a global leader in the female aesthetics industry, GC Aesthetics not only provides products across our Eurosilicone® and Nagor® brands, but also delivers qualified information to help patient's make informed choices.

All dimensions and weights are approximate and slight variations may occur between the value stated in the product catalogue and the actual dimensions and/or weight of the product.



Our Brands

EUROSILICONE®

French manufacturer with 30 years of experience in the implants market.

NAGOR®

The only British manufacturer with 39 years of experience in the implants market.

Silgel[®]

A topical silicone gel designed to minimise the appearance of scars.

GC Aesthetics products intended for sale in Europe are CE certified according to the applicable European medical device regulations.




0086



Our Breast Implants Portfolio at a Glance...


Our Breast Implants Portfolio at a Glance

Round Breast Implants



Please note that this scale is subjective, based on the feel of the product, and is not based on actual test results.

Anatomical Breast Implants



Please note that this scale is subjective, based on the feel of the product, and is not based on actual test results.

Contents

Breast Aesthetics

Round Implants	5
Anatomical Implants	25
Sizers	37
Tissue Expanders	39

Body Aesthetics

Gluteal Implants	50
Calf Implants	54
Testicular Implants	55
Facial Implants	57

Wound Management

Silgel™	62
---------	----

GCA Comfort Guarantee

GCA Comfort Guarantee	66
-----------------------	----

Dimensions & Volumes Chart

Round Collection™	70
IMPLEO™	73

Round Implants

Eurosilicone®

Round Collection™	6
-------------------	---

Nagor®

IMPLEO™	17
---------	----



Round Collection™


by EUROSILICONE®

Product Overview


- ✓ 12 different configurations with 233 variations.
- ✓ **Variable projections:** Low, medium, medium extra-high, high, extra high – which meet the patient desire for low, moderate or large augmentation.
- ✓ **Shell surface:** Smooth or Cristalline™ Micro-textured.
- ✓ **Silicone gel ParaGel™:** Soft Cohesive or Natural Cohesive.
- ✓ **ParagelGard™ 360° Barrier Layer:** Eurosilicone's barrier layer technology minimising gel diffusion.




Natural Cohesive and Micro-textured



Soft Cohesive and Micro-textured



Soft Cohesive and Smooth



Please refer to the pull-out section at the back for dimensions and volumes of all the Eurosilicone® Round Collection™ Implants.

Round Collection™

by EUROSILICONE®


Cristalline Paragel

Natural Cohesive and Microtextured



ES 802N Medium Profile	ES 805N Medium - Extra High Profile	ES 810N High Profile	ES 812N Extra High Profile

DIAMETER / WIDTH	PROJECTION	LVC




ES 802N					Medium Profile	
Product Code	Diameter (cm)	Projection (cm)	Volume (g/cc)	LVC (cm)	Sizer (Gel-filled)	Sizer (Inflatable)
802100N	8.5	2.8	100	5.6	35R100	30100
802120N	9.0	2.9	120	5.9	35R120	30120
802140N	9.3	3.2	140	6.2	35R140	30140
802160N	10.0	3.1	160	6.4	35R160	30160
802180N	10.6	3.2	180	6.8	35R180	30180
802200N	11.1	3.3	200	7.1	35R200	30200
802220N	11.5	3.2	220	7.2	35R220	30220
802240N	11.7	3.5	240	7.5	35R240	30240
802260N	11.9	3.6	260	7.6	35R260	30260
802280N	12.3	3.7	280	7.9	35R280	30280
802300N	12.4	3.9	300	8.0	35R300	30300
802325N	12.7	4.0	325	8.2	35R325	30325
802350N	13.0	4.0	350	8.4	35R350	30350
802375N	13.3	4.2	375	8.6	35R375	30375
802400N	13.7	4.2	400	8.8	35R400	30400
802450N	14.1	4.4	450	9.1	35R450	30450
802500N	14.5	4.6	500	9.4		
802550N	15.0	4.8	550	9.8		
802600N	15.8	4.5	600	9.9		
802650N	15.9	4.9	650	10.2		
802700N	15.9	5.4	700	10.5		
802750N	15.9	5.8	750	10.9		
802800N	15.8	6.0	800	11.0		




ES 805N					Medium - Extra High Profile	
Product Code	Diameter (cm)	Projection (cm)	Volume (g/cc)	LVC (cm)		
805170N	9.9	3.4	170	6.6		
805195N	10.5	3.5	195	6.9		
805215N	10.6	3.7	215	7.1		
805240N	11.4	3.7	240	7.5		
805265N	11.4	3.9	265	7.6		
805290N	11.8	4.0	290	7.8		
805310N	12.0	4.1	310	8.0		
805335N	12.1	4.3	335	8.2		
805360N	12.5	4.4	360	8.4		
805390N	12.8	4.6	390	8.7		
805420N	13.0	4.7	420	8.9		
805450N	13.2	4.8	450	9.0		
805480N	13.9	4.9	480	9.4		

Depending on the product specification, 1g is approximately equivalent to 1cc.

Please note Inflatable Sizers are close but not exact to the measurement of indicated implant.



ES 810N High Profile				
Product Code	Diameter (cm)	Projection (cm)	Volume (g/cc)	LVC (cm)
810155N	9.0	3.6	155	6.4
810180N	9.5	3.8	180	6.8
810215N	10.0	4.0	215	7.1
810240N	10.3	4.2	240	7.4
810265N	10.6	4.4	265	7.7
810295N	10.9	4.6	295	7.9
810320N	11.2	4.8	320	8.2
810350N	11.5	5.0	350	8.5
810390N	12.0	5.2	390	8.8
810435N	12.5	5.4	435	9.2
810500N	13.0	5.8	500	9.7
810575N	13.5	6.2	575	10.2
810665N	14.0	6.6	665	10.7
810765N	14.5	7.0	765	11.2



ES 812N Extra High Profile						
Product Code	Diameter (cm)	Projection (cm)	Volume (g/cc)	LVC (cm)	Sizer (Gel-filled)	Sizer (Inflatable)
812100N	7.3	3.4	100	5.5		33100
812120N	8.0	3.5	120	5.9		33120
812140N	8.5	3.8	140	6.3		33140
812160N	8.8	3.9	160	6.5		
812180N	9.1	4.3	180	7.0	GFS-RCEH 01	33180
812200N	9.4	4.3	200	7.1	GFS-RCEH 02	33200
812220N	9.7	4.4	220	7.3	GFS-RCEH 03	33220
812240N	10.2	4.4	240	7.5	GFS-RCEH 04	33240
812260N	10.4	4.5	260	7.6	GFS-RCEH 05	33260
812280N	10.6	4.7	280	7.9	GFS-RCEH 06	33280
812300N	10.8	4.9	300	8.1	GFS-RCEH 07	33300
812325N	11.0	5.1	325	8.3	GFS-RCEH 08	33325
812350N	11.3	5.1	350	8.4	GFS-RCEH 09	33350
812375N	11.6	5.2	375	8.6	GFS-RCEH 10	33375
812400N	11.8	5.5	400	9.0	GFS-RCEH 11	33400
812450N	12.1	5.7	450	9.2	GFS-RCEH 12	33450
812500N	12.3	6.1	500	9.6		
812550N	12.5	6.4	550	9.9		
812600N	13.1	6.4	600	10.2		
812650N	13.6	6.4	650	10.4		
812700N	13.4	7.1	700	10.8		
812750N	13.5	7.4	750	11.1		
812800N	13.6	7.7	800	11.4		

Depending on the product specification, 1g is approximately equivalent to 1cc.
Please note Inflatable Sizers are close but not exact to the measurement of indicated implant.

Round Collection™

by EUROSILICONE®

Cristalline Paragel
Soft Cohesive and
Microtextured



ES 80 Low Profile	ES 801 Medium Profile	ES 81 High Profile	ES 811 Extra High Profile

 DIAMETER / WIDTH	 PROJECTION	 LVC
----------------------	----------------	---------





ES 80 Low Profile 85% Fill						
Product Code	Diameter (cm)	Projection (cm)	Volume (g/cc)	LVC (cm)	Sizer (Gel-filled)	Sizer (Inflatable)
80100	9.6	2.0	100	5.6	35M100	30100
80120	9.7	2.4	120	5.9	35M120	30120
80140	10.2	2.5	140	6.1	35M140	30140
80160	10.4	2.8	160	6.4	35M160	30160
80180	10.9	2.9	180	6.7	35M180	30180
80200	11.6	2.8	200	7.0	35M200	30200
80220	11.8	3.0	220	7.2	35M220	30220
80240	12.0	3.1	240	7.3	35M240	30240
80260	12.2	3.3	260	7.5	35M260	30260
80280	12.5	3.3	280	7.7	35M280	30280
80300	12.9	3.4	300	7.9	35M300	30300
80325	13.2	3.5	325	8.1	35M325	30325
80350	13.5	3.6	350	8.3	35M350	30350
80375	13.8	3.7	375	8.5	35M375	30375
80400	14.2	3.7	400	8.7	35M400	30400
80450	14.7	3.8	450	9.0	35M450	30450
80500	15.3	4.0	500	9.4		



ES 81 High Profile 85% Fill						
Product Code	Diameter (cm)	Projection (cm)	Volume (g/cc)	LVC (cm)	Sizer (Gel-filled)	Sizer (Inflatable)
81100	7.9	2.9	100	5.4		33100
81120	8.9	2.8	120	5.8		33120
81140	8.8	3.2	140	6.0		33140
81160	9.3	3.3	160	6.3		
81180	9.5	3.6	180	6.6	GFS-RCH 01	33180
81200	9.9	3.7	200	6.8	GFS-RCH 02	33200
81220	10.3	3.7	220	7.0	GFS-RCH 03	33220
81240	10.5	3.9	240	7.2	GFS-RCH 04	33240
81260	10.7	4.1	260	7.5	GFS-RCH 05	33260
81280	10.9	4.3	280	7.7	GFS-RCH 06	33280
81300	11.1	4.4	300	7.8	GFS-RCH 07	33300
81325	11.4	4.5	325	8.0	GFS-RCH 08	33325
81350	11.7	4.7	350	8.3	GFS-RCH 09	33350
81375	11.7	4.9	375	8.5	GFS-RCH 10	33375
81400	12.2	4.9	400	8.7	GFS-RCH 11	33400
81450	12.5	5.3	450	9.1	GFS-RCH 12	33450
81500	12.6	5.8	500	9.5		



ES 801 Medium Profile 100% Fill						
Product Code	Diameter (cm)	Projection (cm)	Volume (g/cc)	LVC (cm)	Sizer (Gel-filled)	Sizer (Inflatable)
801100	8.5	2.8	100	5.6	35R100	30100
801120	9.0	2.9	120	5.9	35R120	30120
801140	9.3	3.2	140	6.2	35R140	30140
801160	10.0	3.1	160	6.4	35R160	30160
801180	10.6	3.2	180	6.8	35R180	30180
801200	11.1	3.3	200	7.1	35R200	30200
801220	11.5	3.2	220	7.2	35R220	30220
801240	11.7	3.5	240	7.5	35R240	30240
801260	11.9	3.6	260	7.6	35R260	30260
801280	12.3	3.7	280	7.9	35R280	30280
801300	12.4	3.9	300	8.0	35R300	30300
801325	12.7	4.0	325	8.2	35R325	30325
801350	13.0	4.0	350	8.4	35R350	30350
801375	13.3	4.2	375	8.6	35R375	30375
801400	13.7	4.2	400	8.8	35R400	30400
801450	14.1	4.4	450	9.1	35R450	30450
801500	14.5	4.6	500	9.4		
801550	15.0	4.8	550	9.8		
801600	15.8	4.5	600	9.9		
801650	15.8	4.9	650	10.2		
801700	15.9	5.4	700	10.5		
801750	15.9	5.8	750	10.9		
801800	15.8	6.0	800	11.0		



ES 811 Extra High Profile 100% Fill						
Product Code	Diameter (cm)	Projection (cm)	Volume (g/cc)	LVC (cm)	Sizer (Gel-filled)	Sizer (Inflatable)
811100	7.3	3.4	100	5.5		33100
811120	8.0	3.5	120	5.9		33120
811140	8.5	3.8	140	6.3		33140
811160	8.8	3.9	160	6.5		
811180	9.1	4.3	180	7.0	GFS-RCEH 01	33180
811200	9.4	4.3	200	7.1	GFS-RCEH 02	33200
811220	9.7	4.4	220	7.3	GFS-RCEH 03	33220
811240	10.2	4.4	240	7.5	GFS-RCEH 04	33240
811260	10.4	4.5	260	7.6	GFS-RCEH 05	33260
811280	10.6	4.7	280	7.9	GFS-RCEH 06	33280
811300	10.8	4.9	300	8.1	GFS-RCEH 07	33300
811325	11.0	5.1	325	8.3	GFS-RCEH 08	33325
811350	11.3	5.1	350	8.4	GFS-RCEH 09	33350
811375	11.6	5.2	375	8.6	GFS-RCEH 10	33375
811400	11.8	5.5	400	9.0	GFS-RCEH 11	33400
811450	12.1	5.7	450	9.2	GFS-RCEH 12	33450
811500	12.3	6.1	500	9.6		
811550	12.5	6.4	550	9.9		
811600	13.1	6.4	600	10.2		
811650	13.6	6.4	650	10.4		
811700	13.4	7.1	700	10.8		
811750	13.5	7.4	750	11.1		
811800	13.6	7.7	800	11.4		

Depending on the product specification, 1g is approximately equivalent to 1cc.
Please note Inflatable Sizers are close but not exact to the measurement of indicated implant.

Depending on the product specification, 1g is approximately equivalent to 1cc.
Please note Inflatable Sizers are close but not exact to the measurement of indicated implant.


Round Collection™

by EUROSILICONE®

Paragel Soft Cohesive and Smooth



ES 26 Low Profile	ES 261 Medium Profile	ES 28 High Profile	ES 281 Extra High Profile



ES 26 Low Profile 85% Fill						
Product Code	Diameter (cm)	Projection (cm)	Volume (g/cc)	LVC (cm)	Sizer (Gel-filled)	Sizer (Inflatable)
26100	9.6	2.0	100	5.6	35M100	30100
26120	9.7	2.4	120	5.9	35M120	30120
26140	10.2	2.5	140	6.1	35M140	30140
26160	10.4	2.8	160	6.4	35M160	30160
26180	10.9	2.9	180	6.7	35M180	30180
26200	11.6	2.8	200	7.0	35M200	30200
26220	11.8	3.0	220	7.2	35M220	30220
26240	12.0	3.1	240	7.3	35M240	30240
26260	12.2	3.3	260	7.5	35M260	30260
26280	12.5	3.3	280	7.7	35M280	30280
26300	12.9	3.4	300	7.9	35M300	30300
26325	13.2	3.5	325	8.1	35M325	30325
26350	13.5	3.6	350	8.3	35M350	30350
26375	13.8	3.7	375	8.5	35M375	30375
26400	14.2	3.7	400	8.7	35M400	30400
26450	14.7	3.8	450	9.0	35M450	30450
26500	15.3	4.0	500	9.4		

Depending on the product specification, 1g is approximately equivalent to 1cc.
Please note Inflatable Sizers are close but not exact to the measurement of indicated implant.




ES 261 Medium Profile 100% Fill						
Product Code	Diameter (cm)	Projection (cm)	Volume (g/cc)	LVC (cm)	Sizer (Gel-filled)	Sizer (Inflatable)
261100	8.5	2.8	100	5.6	35R100	30100
261120	9.0	2.9	120	5.9	35R120	30120
261140	9.3	3.2	140	6.2	35R140	30140
261160	10.0	3.1	160	6.4	35R160	30160
261180	10.6	3.2	180	6.8	35R180	30180
261200	11.1	3.3	200	7.1	35R200	30200
261220	11.5	3.2	220	7.2	35R220	30220
261240	11.7	3.5	240	7.5	35R240	30240
261260	11.9	3.6	260	7.6	35R260	30260
261280	12.3	3.7	280	7.9	35R280	30280
261300	12.4	3.9	300	8.0	35R300	30300
261325	12.7	4.0	325	8.2	35R325	30325
261350	13.0	4.0	350	8.4	35R350	30350
261375	13.3	4.2	375	8.6	35R375	30375
261400	13.7	4.2	400	8.8	35R400	30400
261450	14.1	4.4	450	9.1	35R450	30450
261500	14.5	4.6	500	9.4		
261550	15.0	4.8	550	9.8		
261600	15.8	4.5	600	9.9		
261650	15.8	4.9	650	10.2		
261700	15.8	5.4	700	10.5		
261750	15.9	5.8	750	10.9		
261800	15.8	6.0	800	11.0		



ES 28 High Profile 85% Fill						
Product Code	Diameter (cm)	Projection (cm)	Volume (g/cc)	LVC (cm)	Sizer (Gel-filled)	Sizer (Inflatable)
28100	7.9	2.9	100	5.4		33100
28120	8.9	2.8	120	5.8		33120
28140	8.8	3.2	140	6.0		33140
28160	9.3	3.3	160	6.3		
28180	9.5	3.6	180	6.6	GFS-RCH 01	33180
28200	9.9	3.7	200	6.8	GFS-RCH 02	33200
28220	10.3	3.7	220	7.0	GFS-RCH 03	33220
28240	10.5	3.9	240	7.2	GFS-RCH 04	33240
28260	10.7	4.1	260	7.5	GFS-RCH 05	33260
28280	10.9	4.3	280	7.7	GFS-RCH 06	33280
28300	11.1	4.4	300	7.8	GFS-RCH 07	33300
28325	11.4	4.5	325	8.0	GFS-RCH 08	33325
28350	11.7	4.7	350	8.3	GFS-RCH 09	33350
28375	11.7	4.9	375	8.5	GFS-RCH 10	33375
28400	12.2	4.9	400	8.7	GFS-RCH 11	33400
28450	12.5	5.3	450	9.1	GFS-RCH 12	33450
28500	12.6	5.8	500	9.5		

Depending on the product specification, 1g is approximately equivalent to 1cc.
Please note Inflatable Sizers are close but not exact to the measurement of indicated implant.



ES 281 Extra High Profile 100% Fill						
Product Code	Diameter (cm)	Projection (cm)	Volume (g/cc)	LVC (cm)	Sizer (Gel-filled)	Sizer (Inflatable)
281100	7.3	3.4	100	5.5		33100
281120	8.0	3.5	120	5.9		33120
281140	8.5	3.8	140	6.3		33140
281160	8.8	3.9	160	6.5		
281180	9.1	4.3	180	7.0	GFS-RCEH 01	33180
281200	9.4	4.3	200	7.1	GFS-RCEH 02	33200
281220	9.7	4.4	220	7.3	GFS-RCEH 03	33220
281240	10.2	4.4	240	7.5	GFS-RCEH 04	33240
281260	10.4	4.5	260	7.6	GFS-RCEH 05	33260
281280	10.6	4.7	280	7.9	GFS-RCEH 06	33280
281300	10.8	4.9	300	8.1	GFS-RCEH 07	33300
281325	11.0	5.1	325	8.3	GFS-RCEH 08	33325
281350	11.3	5.1	350	8.4	GFS-RCEH 09	33350
281375	11.6	5.2	375	8.5	GFS-RCEH 10	33375
281400	11.8	5.5	400	9.0	GFS-RCEH 11	33400
281450	12.1	5.7	450	9.2	GFS-RCEH 12	33450
281500	12.3	6.1	500	9.6		
281550	12.5	6.4	550	9.9		
281600	13.1	6.4	600	10.2		
281650	13.6	6.4	650	10.4		
281700	13.4	7.1	700	10.8		
281750	13.5	7.4	750	11.1		
281800	13.6	7.7	800	11.4		

IMPLEO™

by NAGÔR®

Product Overview




- ✓ Smooth or Nagotex™ Mid-textured surface.
- ✓ SiloGard™ 360° Barrier Shell.
- ✓ Extensive choice of fully filled profiles with 38 size options.
- ✓ SiloGel Twist™: 6th generation soft form stable high cohesive silicone gel.




IMPLEO™

by NAGÔR®

Soft Cohesive and Mid-textured




IMP-MR					Moderate Profile	
Product Code	Base Width (cm)	Anterior Proj. (cm)	Volume (g/cc)	LVC (cm)	Sizer (Gel-filled)	Sizer (Inflatable)
IMP-MR 150	9.4	3.4	165	6.5		
IMP-MR 180	10.0	3.5	190	6.8		
IMP-MR 210	10.5	3.6	220	7.1	GS-MR 210	INF-RS 21
IMP-MR 240	11.0	3.9	255	7.4	GS-MR 240	INF-RS 22
IMP-MR 270	11.2	4.1	278	7.7	GS-MR 270	INF-RS 23
IMP-MR 300	11.6	4.2	300	7.9	GS-MR 300	INF-RS 24
IMP-MR 330	12.0	4.3	340	8.1	GS-MR 330	INF-RS 25
IMP-MR 360	12.3	4.4	365	8.3	GS-MR 360	INF-RS 26
IMP-MR 390	12.7	4.4	390	8.5	GS-MR 390	INF-RS 27
IMP-MR 420	13.1	4.6	430	8.8	GS-MR 420	INF-RS 28
IMP-MR 460	13.6	4.7	480	9.1	GS-MR 460	INF-RS 29



IMP-HR					High Profile	
Product Code	Base Width (cm)	Anterior Proj. (cm)	Volume (g/cc)	LVC (cm)	Sizer (Gel-filled)	Sizer (Inflatable)
IMP-HR 150	8.4	4.0	150	6.5		INF-RS 02
IMP-HR 180	8.9	4.2	180	6.8		INF-RS 03
IMP-HR 210	9.5	4.2	210	7.1	GS-HR 210	INF-RS 05
IMP-HR 240	10.0	4.4	240	7.4	GS-HR 240	INF-RS 06
IMP-HR 270	10.2	4.8	270	7.8	GS-HR 270	INF-RS 07
IMP-HR 300	10.7	4.9	300	8.1	GS-HR 300	INF-RS 09
IMP-HR 330	11.2	4.9	330	8.3	GS-HR 330	INF-RS 10
IMP-HR 360	11.6	5.0	360	8.5	GS-HR 360	INF-RS 11
IMP-HR 390	11.8	5.1	375	8.6	GS-HR 390	INF-RS 12
IMP-HR 420	12.4	5.1	425	8.9	GS-HR 420	INF-RS 14
IMP-HR 460	12.7	5.4	460	9.2	GS-HR 460	INF-RS 15
IMP-HR 520	12.9	5.6	505	9.4		INF-RS 16
IMP-HR 560	13.3	5.7	550	9.7		INF-RS 17

Depending on the product specification, 1g is approximately equivalent to 1cc.
Please note Inflatable Sizers are close but not exact to the measurement of indicated implant.



IMP-EHR					Extra High Profile	
Product Code	Base Width (cm)	Anterior Proj. (cm)	Volume (g/cc)	LVC (cm)	Size (Gel-filled)	Size (Inflatable)
IMP-EHR 150	8.2	4.3	161	6.7		
IMP-EHR 180	8.4	4.6	180	7.0		
IMP-EHR 210	9.0	4.8	210	7.3	GS-EHR 210	INF-RS 30
IMP-EHR 240	9.4	4.8	235	7.5	GS-EHR 240	INF-RS 31
IMP-EHR 270	9.8	4.9	260	7.7	GS-EHR 270	INF-RS 32
IMP-EHR 300	10.4	5.0	291	8.1	GS-EHR 300	INF-RS 33
IMP-EHR 330	10.9	5.2	332	8.4	GS-EHR 330	INF-RS 34
IMP-EHR 360	11.2	5.4	360	8.7	GS-EHR 360	INF-RS 35
IMP-EHR 390	11.5	5.5	400	8.9	GS-EHR 390	INF-RS 36
IMP-EHR 420	11.9	5.6	435	9.1	GS-EHR 420	INF-RS 37
IMP-EHR 460	12.2	5.8	460	9.3	GS-EHR 460	INF-RS 38
IMP-EHR 520	12.6	6.1	520	9.7		
IMP-EHR 560	13.0	6.4	575	10.2		
IMP-EHR 620	13.3	6.6	620	10.4		

Depending on the product specification, 1g is approximately equivalent to 1cc.
Please note Inflatable Sizes are close but not exact to the measurement of indicated implant.

Please refer to the pull-out section at the back for dimensions and volumes of all Nagor® IMPLEO™ implants.

IMPLEO Smooth™

by NAGÔR®


Soft Cohesive and Smooth




IMP-SMR Moderate Profile	IMP-SHR High Profile	IMP-SEHR Extra High Profile

 DIAMETER / WIDTH	 PROJECTION	 LVC
----------------------	----------------	---------






IMP-SMR Moderate Profile						
Product Code	Base Width (cm)	Anterior Proj. (cm)	Volume (g/cc)	LVC (cm)	Sizer (Gel-filled)	Sizer (Inflatable)
IMP-SMR 150	9.4	3.4	165	6.5		
IMP-SMR 180	10.0	3.5	190	6.8		
IMP-SMR 210	10.5	3.6	220	7.1	GS-MR 210	INF-RS 21
IMP-SMR 240	11.0	3.9	255	7.4	GS-MR 240	INF-RS 22
IMP-SMR 270	11.2	4.1	278	7.7	GS-MR 270	INF-RS 23
IMP-SMR 300	11.6	4.2	300	7.9	GS-MR 300	INF-RS 24
IMP-SMR 330	12.0	4.3	340	8.1	GS-MR 330	INF-RS 25
IMP-SMR 360	12.3	4.4	365	8.3	GS-MR 360	INF-RS 26
IMP-SMR 390	12.7	4.4	390	8.5	GS-MR 390	INF-RS 27
IMP-SMR 420	13.1	4.6	430	8.8	GS-MR 420	INF-RS 28
IMP-SMR 460	13.6	4.7	480	9.1	GS-MR 460	INF-RS 29



IMP-SHR High Profile						
Product Code	Base Width (cm)	Anterior Proj. (cm)	Volume (g/cc)	LVC (cm)	Sizer (Gel-filled)	Sizer (Inflatable)
IMP-SHR 150	8.4	4.0	150	6.5		INF-RS 02
IMP-SHR 180	8.9	4.2	180	6.8		INF-RS 03
IMP-SHR 210	9.5	4.2	210	7.1	GS-HR 210	INF-RS 05
IMP-SHR 240	10.0	4.4	240	7.4	GS-HR 240	INF-RS 06
IMP-SHR 270	10.2	4.8	270	7.8	GS-HR 270	INF-RS 07
IMP-SHR 300	10.7	4.9	300	8.1	GS-HR 300	INF-RS 09
IMP-SHR 330	11.2	4.9	330	8.3	GS-HR 330	INF-RS 10
IMP-SHR 360	11.6	5.0	360	8.5	GS-HR 360	INF-RS 11
IMP-SHR 390	11.8	5.1	375	8.6	GS-HR 390	INF-RS 12
IMP-SHR 420	12.4	5.1	425	8.9	GS-HR 420	INF-RS 14
IMP-SHR 460	12.7	5.4	460	9.2	GS-HR 460	INF-RS 15
IMP-SHR 520	12.9	5.6	505	9.4		INF-RS 16
IMP-SHR 560	13.3	5.7	550	9.7		INF-RS 17

Depending on the product specification, 1g is approximately equivalent to 1cc.
Please note Inflatable Sizers are close but not exact to the measurement of indicated implant.



IMP-SEHR Extra High Profile						
Product Code	Base Width (cm)	Anterior Proj. (cm)	Volume (g/cc)	LVC (cm)	Sizer (Gel-filled)	Sizer (Inflatable)
IMP-SEHR 150	8.2	4.3	161	6.7		
IMP-SEHR 180	8.4	4.6	180	7.0		
IMP-SEHR 210	9.0	4.8	210	7.3	GS-EHR 210	INF-RS 30
IMP-SEHR 240	9.4	4.8	235	7.5	GS-EHR 240	INF-RS 31
IMP-SEHR 270	9.8	4.9	260	7.7	GS-EHR 270	INF-RS 32
IMP-SEHR 300	10.4	5.0	291	8.1	GS-EHR 300	INF-RS 33
IMP-SEHR 330	10.9	5.2	332	8.4	GS-EHR 330	INF-RS 34
IMP-SEHR 360	11.2	5.4	360	8.7	GS-EHR 360	INF-RS 35
IMP-SEHR 390	11.5	5.5	400	8.9	GS-EHR 390	INF-RS 36
IMP-SEHR 420	11.9	5.6	435	9.1	GS-EHR 420	INF-RS 37
IMP-SEHR 460	12.2	5.8	460	9.3	GS-EHR 460	INF-RS 38
IMP-SEHR 520	12.6	6.1	520	9.7		
IMP-SEHR 560	13.0	6.4	575	10.2		
IMP-SEHR 620	13.3	6.6	620	10.4		

Depending on the product specification, 1g is approximately equivalent to 1cc.
Inflatable Sizers available on page 42.

Anatomical Implants

Eurosilicone®

The Matrix™	26
-------------	----

Nagor®

CoGel®	31
--------	----



The Matrix™

by EUROSILICONE®

Designed for
Performance

Natural Cohesive and Micro-textured

Product Overview

- ✓ Tapered edges for less visibility and palpability.
- ✓ Firm upper pole for shape maintenance.
- ✓ Higher projection for added lift.
- ✓ ParaGel™ form stable high cohesive silicone gel.
- ✓ Crystalline Micro-textured.
- ✓ ParagelGard™ 360° Barrier Layer.
- ✓ 7 orientation marker dots.
- ✓ Gel-filled and Inflatable
Sizes available for core sizes.



TM F1 (F) Full Height (1) Low Projection						
Product Code	Width (cm)	Height (cm)	Proj. (cm)	Arc Lgth. (cm)	Volume (cc)	LVC (cm)
TM F1-120	9.5	10.0	3.0	4.6	120	6.9
TM F1-140	10.0	10.5	3.2	4.8	140	7.3
TM F1-170	10.5	11.0	3.4	5.1	170	7.7
TM F1-195	11.0	11.5	3.6	5.4	195	8.1
TM F1-225	11.5	12.0	3.8	5.6	225	8.5
TM F1-250	12.0	12.5	4.0	5.9	250	8.9
TM F1-290	12.5	13.0	4.2	6.2	290	9.3
TM F1-320	13.0	13.5	4.4	6.4	320	9.7
TM F1-370	13.5	14.0	4.6	6.7	370	10.1
TM F1-420	14.0	14.5	4.8	7.0	420	10.5

TM M1 (M) Moderate Height (1) Low Projection						
Product Code	Width (cm)	Height (cm)	Proj. (cm)	Arc Lgth. (cm)	Volume (cc)	LVC (cm)
TM M1-110	9.5	9.0	3.0	4.4	110	6.5
TM M1-125	10.0	9.5	3.2	4.6	125	6.9
TM M1-145	10.5	10.0	3.4	4.9	145	7.3
TM M1-175	11.0	10.5	3.6	5.2	175	7.7
TM M1-195	11.5	11.0	3.8	5.4	195	8.1
TM M1-225	12.0	11.5	4.0	5.7	225	8.5
TM M1-255	12.5	12.0	4.2	6.0	255	8.9
TM M1-295	13.0	12.5	4.4	6.2	295	9.3
TM M1-330	13.5	13.0	4.6	6.5	330	9.7
TM M1-375	14.0	13.5	4.8	6.8	375	10.1

TM L1 (L) Low Height (1) Low Projection						
Product Code	Width (cm)	Height (cm)	Proj. (cm)	Arc Lgth. (cm)	Volume (cc)	LVC (cm)
TM L1-100	9.5	8.0	3.0	4.2	100	6.1
TM L1-115	10.0	8.5	3.2	4.4	115	6.5
TM L1-130	10.5	9.0	3.4	4.7	130	6.9
TM L1-160	11.0	9.5	3.6	5.0	160	7.3
TM L1-180	11.5	10.0	3.8	5.2	180	7.7
TM L1-210	12.0	10.5	4.0	5.5	210	8.1
TM L1-245	12.5	11.0	4.2	5.8	245	8.5
TM L1-280	13.0	11.5	4.4	6.1	280	8.9
TM L1-320	13.5	12.0	4.6	6.3	320	9.3
TM L1-355	14.0	12.5	4.8	6.6	355	9.7

TM F2 (F) Full Height (2) Moderate Projection								
Product Code	Width (cm)	Height (cm)	Proj. (cm)	Arc Lgth. (cm)	Volume (cc)	LVC (cm)	Sizer (Gel-filled)	Sizer (Inflatable)
TM F2-150	9.5	10.0	3.7	5.2	150	7.6	GFS-TMF2 01	INF-TMF2S 01
TM F2-175	10.0	10.5	3.9	5.4	175	8.0	GFS-TMF2 02	INF-TMF2S 02
TM F2-195	10.5	11.0	4.1	5.7	195	8.4	GFS-TMF2 03	INF-TMF2S 03
TM F2-225	11.0	11.5	4.3	6.0	225	8.8	GFS-TMF2 04	INF-TMF2S 04
TM F2-250	11.5	12.0	4.5	6.2	250	9.2	GFS-TMF2 05	INF-TMF2S 05
TM F2-285	12.0	12.5	4.7	6.5	285	9.6	GFS-TMF2 06	INF-TMF2S 06
TM F2-320	12.5	13.0	4.9	6.8	320	10.0	GFS-TMF2 07	INF-TMF2S 07
TM F2-360	13.0	13.5	5.1	7.1	360	10.4	GFS-TMF2 08	INF-TMF2S 08
TM F2-410	13.5	14.0	5.3	7.3	410	10.8	GFS-TMF2 09	INF-TMF2S 09
TM F2-465	14.0	14.5	5.5	7.6	465	11.2	GFS-TMF2 10	INF-TMF2S 10
TM F2-540	14.5	15.0	5.7	7.9	540	11.6		
TM F2-615	15.0	15.5	5.9	8.1	615	12.0		
TM F2-670	15.5	16.0	6.1	8.4	670	12.4		

TM M2 (M) Moderate Height (2) Moderate Projection								
Product Code	Width (cm)	Height (cm)	Proj. (cm)	Arc Lgth. (cm)	Volume (cc)	LVC (cm)	Sizer (Gel-filled)	Sizer (Inflatable)
TM M2-140	9.5	9.0	3.7	5.0	140	7.3	GFS-TM M2 01	INF-TMM2S 01
TM M2-160	10.0	9.5	3.9	5.2	160	7.6	GFS-TMM2 02	INF-TMM2S 02
TM M2-180	10.5	10.0	4.1	5.5	180	8.0	GFS-TMM2 03	INF-TMM2S 03
TM M2-210	11.0	10.5	4.3	5.8	210	8.4	GFS-TMM2 04	INF-TMM2S 04
TM M2-235	11.5	11.0	4.5	6.1	235	8.8	GFS-TMM2 05	INF-TMM2S 05
TM M2-270	12.0	11.5	4.7	6.3	270	9.2	GFS-TMM2 06	INF-TMM2S 06
TM M2-305	12.5	12.0	4.9	6.6	305	9.6	GFS-TMM2 07	INF-TMM2S 07
TM M2-350	13.0	12.5	5.1	6.9	350	10.0	GFS-TMM2 08	INF-TMM2S 08
TM M2-385	13.5	13.0	5.3	7.1	385	10.4	GFS-TMM2 09	INF-TMM2S 09
TM M2-435	14.0	13.5	5.5	7.4	435	10.8	GFS-TMM2 10	INF-TMM2S 10
TM M2-495	14.5	14.0	5.7	7.7	495	11.2		
TM M2-560	15.0	14.5	5.9	8.0	560	11.6		
TM M2-615	15.5	15.0	6.1	8.2	615	12.0		

TM L2 (L) Low Height (2) Moderate Projection							
Product Code	Width (cm)	Height (cm)	Proj. (cm)	Arc Lgth. (cm)	Volume (cc)	LVC (cm)	
TM L2-120	9.5	8.0	3.7	4.8	120	6.9	
TM L2-145	10.0	8.5	3.9	5.1	145	7.3	
TM L2-175	10.5	9.0	4.1	5.3	175	7.7	
TM L2-195	11.0	9.5	4.3	5.6	195	8.1	
TM L2-225	11.5	10.0	4.5	5.9	225	8.5	
TM L2-260	12.0	10.5	4.7	6.2	260	8.9	
TM L2-295	12.5	11.0	4.9	6.4	295	9.3	
TM L2-330	13.0	11.5	5.1	6.7	330	9.7	
TM L2-375	13.5	12.0	5.3	7.0	375	10.1	
TM L2-415	14.0	12.5	5.5	7.2	415	10.5	

Please note Inflatable Sizers are close but not exact to the measurement of indicated implant.

TM F3 (F) Full Height (3) High Projection							
Product Code	Width (cm)	Height (cm)	Proj. (cm)	Arc Lgth. (cm)	Volume (cc)	LVC (cm)	Sizer (Gel-filled)
TM F3-175	9.5	10.0	4.5	5.9	175	8.5	GFS-TMF3 01
TM F3-205	10.0	10.5	4.7	6.2	205	8.9	GFS-TMF3 02
TM F3-230	10.5	11.0	4.9	6.4	230	9.3	GFS-TMF3 03
TM F3-275	11.0	11.5	5.2	6.8	275	9.8	GFS-TMF3 04
TM F3-305	11.5	12.0	5.4	7.1	305	10.2	GFS-TMF3 05
TM F3-345	12.0	12.5	5.6	7.3	345	10.6	GFS-TMF3 06
TM F3-385	12.5	13.0	5.9	7.7	385	11.1	GFS-TMF3 07
TM F3-430	13.0	13.5	6.2	8.1	430	11.6	GFS-TMF3 08
TM F3-475	13.5	14.0	6.3	8.2	475	11.9	GFS-TMF3 09
TM F3-545	14.0	14.5	6.4	8.4	545	12.2	GFS-TMF3 10
TM F3-610	14.5	15.0	6.6	8.7	610	12.6	
TM F3-680	15.0	15.5	6.8	9.0	680	13.0	
TM F3-730	15.5	16.0	7.0	9.2	730	13.3	

TM M3 (M) Moderate Height (3) High Projection							
Product Code	Width (cm)	Height (cm)	Proj. (cm)	Arc Lgth. (cm)	Volume (cc)	LVC (cm)	Sizer (Gel-filled)
TM M3-160	9.5	9.0	4.5	5.7	160	8.1	GFS-TMM3 01
TM M3-190	10.0	9.5	4.7	6.0	190	8.5	GFS-TMM3 02
TM M3-215	10.5	10.0	4.9	6.3	215	8.9	GFS-TMM3 03
TM M3-255	11.0	10.5	5.2	6.6	255	9.4	GFS-TMM3 04
TM M3-290	11.5	11.0	5.4	6.9	290	9.8	GFS-TMM3 05
TM M3-320	12.0	11.5	5.6	7.2	320	10.2	GFS-TMM3 06
TM M3-370	12.5	12.0	5.9	7.5	370	10.7	GFS-TMM3 07
TM M3-410	13.0	12.5	6.2	7.9	410	11.2	GFS-TMM3 08
TM M3-445	13.5	13.0	6.3	8.1	445	11.5	GFS-TMM3 09
TM M3-515	14.0	13.5	6.4	8.3	515	11.8	GFS-TMM3 10
TM M3-575	14.5	14.0	6.6	8.5	575	12.2	
TM M3-645	15.0	14.5	6.8	8.8	645	12.6	
TM M3-695	15.5	15.0	7.0	9.1	695	13.0	

TM L3 (L) Low Height (1) Low Projection							
Product Code	Width (cm)	Height (cm)	Proj. (cm)	Arc Lgth. (cm)	Volume (cc)	LVC (cm)	
TM L3-145	9.5	8.0	4.5	5.6	145	7.8	
TM L3-175	10.0	8.5	4.7	5.9	175	8.2	
TM L3-200	10.5	9.0	4.9	6.1	200	8.6	
TM L3-240	11.0	9.5	5.2	6.5	240	9.1	
TM L3-270	11.5	10.0	5.4	6.8	270	9.5	
TM L3-310	12.0	10.5	5.6	7.0	310	9.9	
TM L3-360	12.5	11.0	5.9	7.4	360	10.4	
TM L3-390	13.0	11.5	6.2	7.8	390	10.9	
TM L3-450	13.5	12.0	6.3	7.9	450	11.2	
TM L3-495	14.0	12.5	6.4	8.1	495	11.5	

CoGel®

by NAGÔR®


Designed for
Performance

Firm Cohesive and Mid-textured

Product Overview

- ✓ Nagotex® Mid-textured Surface.
- ✓ Tapered edges for less visibility and palpability.
- ✓ Firm upper pole for shape maintenance.
- ✓ Higher projection for added lift.
- ✓ SiloGard™ 360° Barrier.
- ✓ SiloGel™ form stable high cohesive gel.
- ✓ 7 Orientation marker dots.
- ✓ Gel-filled and Inflatable Sizers available for core sizes.





XF1 (F) Full Height (1) Low Projection						
Product Code	Base Width (cm)	Height (cm)	Anterior Proj. (cm)	Arc Lgth. (cm)	Volume (g)	LVC (cm)
XF1-120	9.5	10.0	3.0	4.3	120	7.1
XF1-145	10.0	10.5	3.2	4.6	145	7.5
XF1-170	10.5	11.0	3.4	4.9	170	7.9
XF1-195	11.0	11.5	3.6	5.1	195	8.3
XF1-235	11.5	12.0	3.8	5.4	235	8.7
XF1-275	12.0	12.5	4.0	5.7	275	9.1
XF1-300	12.5	13.0	4.2	5.9	300	9.5
XF1-350	13.0	13.5	4.4	6.2	350	9.9
XF1-390	13.5	14.0	4.6	6.4	390	10.3
XF1-430	14.0	14.5	4.8	6.7	430	10.7

XM1 (M) Moderate Height (1) Low Projection						
Product Code	Base Width (cm)	Height (cm)	Anterior Proj. (cm)	Arc Lgth. (cm)	Volume (g)	LVC (cm)
XM1-105	9.5	9.0	3.0	4.3	105	6.6
XM1-135	10.0	9.5	3.2	4.6	135	7.0
XM1-155	10.5	10.0	3.4	4.9	155	7.4
XM1-185	11.0	10.5	3.6	5.1	185	7.8
XM1-215	11.5	11.0	3.8	5.4	215	8.2
XM1-250	12.0	11.5	4.0	5.7	250	8.6
XM1-275	12.5	12.0	4.2	5.9	275	9.0
XM1-325	13.0	12.5	4.4	6.2	325	9.4
XM1-360	13.5	13.0	4.6	6.4	360	9.8
XM1-410	14.0	13.5	4.8	6.7	410	10.2

XL1 (L) Low Height (1) Low Projection						
Product Code	Base Width (cm)	Height (cm)	Anterior Proj. (cm)	Arc Lgth. (cm)	Volume (g)	LVC (cm)
XL1-105	9.5	8.0	3.0	4.3	105	6.1
XL1-120	10.0	8.5	3.2	4.6	120	6.5
XL1-140	10.5	9.0	3.4	4.9	140	6.9
XL1-170	11.0	9.5	3.6	5.1	170	7.2
XL1-210	11.5	10.0	3.8	5.4	210	7.6
XL1-230	12.0	10.5	4.0	5.7	230	8.0
XL1-250	12.5	11.0	4.2	5.9	250	8.4
XL1-300	13.0	11.5	4.4	6.2	300	8.8
XL1-330	13.5	12.0	4.6	6.4	330	9.1
XL1-380	14.0	12.5	4.8	6.7	380	9.6

XF2 (F) Full Height (2) Moderate Projection								
Product Code	Base Width (cm)	Height (cm)	Anterior Proj. (cm)	Arc Lgth. (cm)	Volume (g)	LVC (cm)	Sizer (Gel-filled)	Sizer (Inflatable)
XF2-145	9.5	10.0	3.7	4.9	145	7.8	GS-XF2 145	INF-XF2S 01
XF2-170	10.0	10.5	3.9	5.2	170	8.1	GS-XF2 170	INF-XF2S 02
XF2-200	10.5	11.0	4.1	5.4	200	8.5	GS-XF2 200	INF-XF2S 03
XF2-230	11.0	11.5	4.3	5.7	230	8.9	GS-XF2 230	INF-XF2S 04
XF2-270	11.5	12.0	4.5	6.0	270	9.3	GS-XF2 270	INF-XF2S 05
XF2-315	12.0	12.5	4.7	6.3	315	9.7	GS-XF2 315	INF-XF2S 06
XF2-350	12.5	13.0	4.9	6.5	350	10.1	GS-XF2 350	INF-XF2S 07
XF2-400	13.0	13.5	5.1	6.8	400	10.5	GS-XF2 400	INF-XF2S 08
XF2-450	13.5	14.0	5.3	7.0	450	10.9	GS-XF2 450	INF-XF2S 09
XF2-505	14.0	14.5	5.5	7.3	505	11.3	GS-XF2 505	INF-XF2S 10

XM2 (M) Moderate Height (2) Moderate Projection								
Product Code	Base Width (cm)	Height (cm)	Anterior Proj. (cm)	Arc Lgth. (cm)	Volume (g)	LVC (cm)	Sizer (Gel-filled)	Sizer (Inflatable)
XM2-130	9.5	9.0	3.7	4.9	130	7.3	GS-XM2 130	INF-XM2S 01
XM2-165	10.0	9.5	3.9	5.2	165	7.6	GS-XM2 165	INF-XM2S 02
XM2-185	10.5	10.0	4.1	5.4	185	8.0	GS-XM2 185	INF-XM2S 03
XM2-215	11.0	10.5	4.3	5.7	215	8.4	GS-XM2 215	INF-XM2S 04
XM2-250	11.5	11.0	4.5	6.0	250	8.8	GS-XM2 250	INF-XM2S 05
XM2-290	12.0	11.5	4.7	6.3	290	9.2	GS-XM2 290	INF-XM2S 06
XM2-320	12.5	12.0	4.9	6.5	320	9.6	GS-XM2 320	INF-XM2S 07
XM2-370	13.0	12.5	5.1	6.8	370	10.0	GS-XM2 370	INF-XM2S 08
XM2-415	13.5	13.0	5.3	7.0	415	10.4	GS-XM2 415	INF-XM2S 09
XM2-445	14.0	13.5	5.5	7.3	445	10.8	GS-XM2 445	INF-XM2S 10

XL2 (L) Low Height (2) Moderate Projection						
Product Code	Base Width (cm)	Height (cm)	Anterior Proj. (cm)	Arc Lgth. (cm)	Volume (g)	LVC (cm)
XL2-120	9.5	8.0	3.7	4.9	120	6.7
XL2-145	10.0	8.5	3.9	5.2	145	7.1
XL2-170	10.5	9.0	4.1	5.4	170	7.5
XL2-195	11.0	9.5	4.3	5.7	195	7.9
XL2-230	11.5	10.0	4.5	6.0	230	8.3
XL2-265	12.0	10.5	4.7	6.3	265	8.7
XL2-295	12.5	11.0	4.9	6.5	295	9.1
XL2-340	13.0	11.5	5.1	6.8	340	9.5
XL2-385	13.5	12.0	5.3	7.0	385	9.9
XL2-425	14.0	12.5	5.5	7.3	425	10.3

XF3 (F) Full Height (3) High Projection								
Product Code	Base Width (cm)	Height (cm)	Anterior Proj. (cm)	Arc Lgth. (cm)	Volume (g)	LVC (cm)	Sizer (Gel-filled)	
XF3-170	9.5	10.0	4.5	5.4	170	8.5	GS-XF3 170	
XF3-200	10.0	10.5	4.7	5.7	200	8.9	GS-XF3 200	
XF3-235	10.5	11.0	4.9	5.9	235	9.3	GS-XF3 235	
XF3-280	11.0	11.5	5.2	6.2	280	9.8	GS-XF3 280	
XF3-320	11.5	12.0	5.4	6.5	320	10.3	GS-XF3 320	
XF3-365	12.0	12.5	5.6	6.8	365	10.7	GS-XF3 365	
XF3-405	12.5	13.0	5.9	7.0	405	11.1	GS-XF3 405	
XF3-465	13.0	13.5	6.2	7.3	465	11.5	GS-XF3 465	
XF3-515	13.5	14.0	6.3	7.5	515	11.8	GS-XF3 515	
XF3-540	14.0	14.5	6.4	7.8	540	12.2	GS-XF3 540	

XM3 (M) Moderate Height (3) High Projection								
Product Code	Base Width (cm)	Height (cm)	Anterior Proj. (cm)	Arc Lgth. (cm)	Volume (g)	LVC (cm)	Sizer (Gel-filled)	
XM3-155	9.5	9.0	4.5	5.4	155	8.1	GS-XM3 155	
XM3-190	10.0	9.5	4.7	5.7	190	8.6	GS-XM3 190	
XM3-215	10.5	10.0	4.9	5.9	215	9.0	GS-XM3 215	
XM3-260	11.0	10.5	5.2	6.2	260	9.4	GS-XM3 260	
XM3-300	11.5	11.0	5.4	6.5	300	9.8	GS-XM3 300	
XM3-335	12.0	11.5	5.6	6.8	335	10.2	GS-XM3 335	
XM3-375	12.5	12.0	5.9	7.0	375	10.6	GS-XM3 375	
XM3-440	13.0	12.5	6.2	7.3	440	11.0	GS-XM3 440	
XM3-485	13.5	13.0	6.3	7.5	485	11.4	GS-XM3 485	
XM3-520	14.0	13.5	6.4	7.8	520	11.8	GS-XM3 520	

XL3 (L) Low Height (1) Low Projection						
Product Code	Base Width (cm)	Height (cm)	Anterior Proj. (cm)	Arc Lgth. (cm)	Volume (g)	LVC (cm)
XL3-160	9.5	8.0	4.5	5.4	160	7.5
XL3-175	10.0	8.5	4.7	5.7	175	8.0
XL3-200	10.5	9.0	4.9	5.9	200	8.4
XL3-240	11.0	9.5	5.2	6.2	240	8.8
XL3-275	11.5	10.0	5.4	6.5	275	9.2
XL3-305	12.0	10.5	5.6	6.8	305	9.6
XL3-350	12.5	11.0	5.9	7.0	350	10.0
XL3-410	13.0	11.5	6.2	7.3	410	10.4
XL3-450	13.5	12.0	6.3	7.5	450	10.8
XL3-480	14.0	12.5	6.4	7.8	480	11.2

Please note Inflatable Sizers are close but not exact to the measurement of indicated implant.

Sizers

Sizers:

General information	38
---------------------	----

Sizers

GCA Aesthetics offers sizers which are indicated to assist with intra-operative volume determination during breast augmentation and reconstruction procedures.

Range includes gel-filled mammary sizers pre-filled with silicone gel or inflatable mammary sizers filled intraoperatively with sterile saline solution.

- ✓ Sterile, single use
- ✓ Available for anatomical and round implants
- ✓ Smooth surface
- ✓ Temporary device – must not be used as a permanent implant
- ✓ For more information on each product's sizers, please refer to that implant's individual product page.

Tissue Expanders

Eurosilicone® and Nagor®

Integrated Valve (IVTE)	44
Remote Valve (RVTE)	48

Integrated Valve Tissue Expanders

by EUROSILICONE® & NAGÔR®



Product Overview

- ✓ 3 different configurations with 16 different variations.
- ✓ Anatomical shape, mid-textured.
- ✓ Re-usable and ergonomically shaped external location device.
- ✓ Large integrated injection site with self-sealing technology which helps to facilitate the transfer of saline solution from the needle to the product.
- ✓ Reinforced base layer to ensure unidirectional expansion.
- ✓ Compliments the GCA anatomical breast implants range (The Matrix™ & CoGel®).

All device dimensions correspond to the suggested volume.

GCA™ Integrated-Valve Tissue Expanders must be filled with a Sterile isotonic saline solution only.


A 21 gauge standard bevel hypodermic needle must be used in order to inflate the device through the injection site.

External Location Device


Supplied sterile with the Integrated-Valve Tissue Expander. Extra non-sterile ELD's can be ordered separately.




EUROSILICONE®



TE-F3 (F) Full Height (3) High Projection				
Product Code	Width (cm)	Height (cm)	Projection (cm)	Volume (cc)
TE-F3 300	11.0	11.5	5.2	300
TE-F3 400	12.0	12.5	5.6	400
TE-F3 500	13.0	13.5	6.2	500
TE-F3 600	14.0	14.5	6.4	600
TE-F3 700	15.0	15.5	6.8	700
TE-F3 800	15.5	16.0	7.0	800




TE-M3 (M) Moderate Height (3) High Projection				
Product Code	Width (cm)	Height (cm)	Projection (cm)	Volume (cc)
TE-M3 200	10.0	9.5	4.7	200
TE-M3 300	11.0	10.5	5.2	300
TE-M3 400	12.0	11.5	5.6	400
TE-M3 500	13.0	12.5	6.2	500
TE-M3 600	14.0	13.5	6.4	600
TE-M3 700	15.0	14.5	6.8	700




TE-L3 (L) Low Height (3) High Projection				
Product Code	Width (cm)	Height (cm)	Projection (cm)	Volume (cc)
TE-L3 200	10.0	8.5	4.7	200
TE-L3 300	11.5	10.0	5.4	300
TE-L3 400	12.5	11.0	5.9	400
TE-L3 500	13.5	12.0	6.3	500


NAGÔR®



TE-XF3 (F) Full Height (3) High Projection				
Product Code	Width (cm)	Height (cm)	Projection (cm)	Volume (ml)
TE-XF3 300	11.0	11.5	5.2	300
TE-XF3 400	12.0	12.5	5.6	400
TE-XF3 500	13.0	13.5	6.2	500
TE-XF3 600	14.0	14.5	6.4	600
TE-XF3 700	15.0	15.5	6.8	700
TE-XF3 800	15.5	16.0	7.0	800



TE-XM3 (M) Moderate Height (3) High Projection				
Product Code	Width (cm)	Height (cm)	Projection (cm)	Volume (ml)
TE-XM3 200	10.0	9.5	4.7	200
TE-XM3 300	11.0	10.5	5.2	300
TE-XM3 400	12.0	11.5	5.6	400
TE-XM3 500	13.0	12.5	6.2	500
TE-XM3 600	14.0	13.5	6.4	600
TE-XM3 700	15.0	14.5	6.8	700



TE-XL3 (L) Low Height (3) High Projection				
Product Code	Width (cm)	Height (cm)	Projection (cm)	Volume (ml)
TE-XL3 200	10.0	8.5	4.7	200
TE-XL3 300	11.5	10.0	5.4	300
TE-XL3 400	12.5	11.0	5.9	400
TE-XL3 500	13.5	12.0	6.3	500

Remote Valve Tissue Expanders

by EUROSILICONE® & NAGÔR®

Available in 5 different smooth configurations:

- ✓ Round shape.
- ✓ Rectangular shape.
- ✓ Cylindrical straight.
- ✓ Cylindrical curved.
- ✓ Croissant shape.

All device dimensions correspond to the suggested fill volume.

GCA™ Remote-Valve Tissue Expanders must be filled with a sterile isotonic solution only.


A fill tube is connected to every tissue expander to allow inflation of the device through an injection port.

All GCA™ Remote-Valve Tissue Expanders with a recommended fill volume of more than 150cc are provided with a 28mm injection port, smaller volume devices are provided with a 20mm injection port.

Round EUROSILICONE®

ES 12

Product Code	Diameter (cm)	Height (cm)	Valve (cm)	Volume (cc)
12010	3.1	2.1	2.0	10
12120	7.1	4.8	2.0	120
12150	9.5	4.0	2.0	150
12200	7.9	5.1	2.8	200
12250	9.0	5.7	2.8	250
12300	10.0	6.3	2.8	300
12375	13.0	5.0	2.8	375
12400	11.0	6.0	2.8	400
12500	12.8	6.3	2.8	500
12600	13.8	6.7	2.8	600
12700	14.8	7.5	2.8	700
12800	15.7	6.7	2.8	800
12900	16.6	6.5	2.8	900
121100	17.4	8.2	2.8	1100




Round NAGÔR®

TE-RD

Product Code	Diameter (cm)	Height (cm)	Valve (cm)	Volume (ml)
TE-RD 010	3.1	2.1	2.0	10
TE-RD 120	7.1	4.8	2.0	120
TE-RD 150	9.5	4.0	2.0	150
TE-RD 200	7.9	5.1	2.8	200
TE-RD 250	9.0	5.7	2.8	250
TE-RD 300	10.0	6.3	2.8	300
TE-RD 375	13.0	5.0	2.8	375
TE-RD 400	11.0	6.0	2.8	400
TE-RD 500	12.8	6.3	2.8	500
TE-RD 600	13.8	6.7	2.8	600
TE-RD 700	14.8	7.5	2.8	700
TE-RD 800	15.7	6.7	2.8	800
TE-RD 900	16.6	6.5	2.8	900
TE-RD 1100	17.4	8.2	2.8	1100

Rectangular EUROSILICONE®

ES 14					
Product Code	Length (cm)	Width (cm)	Projection (cm)	Valve (cm)	Volume (cc)
14016	4.0	2.1	2.1	2.0	16
14035	4.0	2.6	3.6	2.0	35
14050	5.0	3.0	3.6	2.0	50
14100	4.9	3.9	5.1	2.0	100
14150	5.8	4.8	5.0	2.0	150
14200	8.5	5.0	4.9	2.8	200
14300	10.0	6.0	6.3	2.8	300
14400	10.8	5.8	5.9	2.8	400
14520	11.7	7.8	5.4	2.8	520
14650	11.7	9.0	5.8	2.8	650
14720	14.5	7.8	5.8	2.8	720
14750	13.5	8.7	5.9	2.8	750
14800	15.5	8.2	5.8	2.8	800
14950	17.5	8.8	5.5	2.8	950




Rectangular NAGÔR®

TE-RC					
Product Code	Length (cm)	Width (cm)	Projection (cm)	Valve (cm)	Volume (ml)
TE-RC 016	4.0	2.1	2.1	2.0	16
TE-RC 035	4.0	2.6	3.6	2.0	35
TE-RC 050	5.0	3.0	3.6	2.0	50
TE-RC 100	4.9	3.9	5.1	2.0	100
TE-RC 150	5.8	4.8	5.0	2.0	150
TE-RC 200	8.5	5.0	4.9	2.8	200
TE-RC 300	10.0	6.0	6.3	2.8	300
TE-RC 400	10.8	5.8	5.9	2.8	400
TE-RC 520	11.7	7.8	5.4	2.8	520
TE-RC 650	11.7	9.0	5.8	2.8	650
TE-RC 720	14.5	7.8	5.8	2.8	720
TE-RC 750	13.5	8.7	5.9	2.8	750
TE-RC 800	15.5	8.2	5.8	2.8	800
TE-RC 950	17.5	8.8	5.5	2.8	950


Cylindrical Straight EUROSILICONE®

ES 16					
Product Code	Length (cm)	Width (cm)	Projection (cm)	Valve (cm)	Volume (cc)
16060	10.0	3.9	2.7	2.0	60
16100	11.0	3.5	3.5	2.0	100
16150	10.0	5.4	4.0	2.0	150
16200	14.0	4.5	4.5	2.8	200
16250	10.2	7.2	5.4	2.8	250
16300	16.0	5.0	5.0	2.8	300
16400	17.5	6.0	6.0	2.8	400



Cylindrical Straight NAGÔR®

TE-CYS					
Product Code	Length (cm)	Width (cm)	Projection (cm)	Valve (cm)	Volume (ml)
TE-CYS 060	10.0	3.9	2.7	2.0	60
TE-CYS 100	11.0	3.5	3.5	2.0	100
TE-CYS 150	10.0	5.4	4.0	2.0	150
TE-CYS 200	14.0	4.5	4.5	2.8	200
TE-CYS 250	10.2	7.2	5.4	2.8	250
TE-CYS 300	16.0	5.0	5.0	2.8	300
TE-CYS 400	17.5	6.0	6.0	2.8	400



Cylindrical Curved EUROSILICONE®

ES 16C

Product Code	Length (cm)	Width (cm)	Projection (cm)	Valve (cm)	Volume (cc)
16C250	14.5	7.7	4.5	2.8	250

Cylindrical Curved NAGÔR®

TE-CYC

Product Code	Length (cm)	Width (cm)	Projection (cm)	Valve (cm)	Volume (ml)
TE-CYC 250	14.5	7.7	4.5	2.8	250

Croissant EUROSILICONE®

ES 19

Product Code	Length (cm)	Width (cm)	Projection (cm)	Valve (cm)	Volume (cc)
19060	8.5	6.0	3.9	2.0	60
19100	10.8	8.5	5.5	2.0	100
19125	10.0	7.0	4.8	2.0	125
19250	13.0	10.0	5.7	2.8	250
19500	17.0	12.0	7.8	2.8	500

Croissant NAGÔR®

TE-CS

Product Code	Length (cm)	Width (cm)	Projection (cm)	Valve (cm)	Volume (ml)
TE-CS 060	8.5	6.0	3.9	2.0	60
TE-CS 100	10.8	8.5	5.5	2.0	100
TE-CS 125	10.0	7.0	4.8	2.0	125
TE-CS 250	13.0	10.0	5.7	2.8	250
TE-CS 500	17.0	12.0	7.8	2.8	500


Body Aesthetics

Gluteal Implants	50
Calf Implants	54
Testicular Implants	55
Facial Implants	57

Gluteal Implants

by EUROSILICONE®


- ✓ 2 different configurations with 17 different variations.
- ✓ Round or anatomical shape.
- ✓ High Cohesive Gel.
- ✓ Smooth shell surface.



ES 62 Smooth Surface, Round			
Product Code	Diameter (cm)	Projection (cm)	Volume (cc)
62140	8.4	4.1	140
62180	9.1	4.3	180
62200	9.4	4.4	200
62220	9.7	4.5	220
62240	10.3	4.7	240
62260	10.4	4.7	260
62280	10.6	4.9	280
62300	10.8	5.0	300



* Not available in Europe.



ES 60 Smooth Surface, Anatomical				
Product Code	Height (cm)	Width (cm)	Projection (cm)	Volume (cc)
60150	13.1	9.6	2.3	150
60210	14.6	10.7	2.6	210
60230	12.0	11.2	3.5	230
60275	12.4	11.4	3.8	275
60300	13.9	10.9	3.6	300
60325	14.7	11.3	3.7	325
60375	14.7	11.7	4.1	375
60425	15.0	12.4	4.5	425
60480	15.5	13.2	4.8	480




* Not available in Europe.

Gluteal Implants

CURVEO™ by NAGÔR®


- ✓ A highly cohesive, form stable gel for a natural look.
- ✓ Orientation markers allowing for more precision during placement.
- ✓ Optimised dimensional ratio between projection and base width, for enhanced aesthetic results.
- ✓ Orientation markers allowing for more precision.
- ✓ Optimised dimensional ratio between projection and base width



GLT-R Round Gluteal Implant			
Product Code	Diameter (mm) ±2	Projection (mm) ±2	Volume (g) ±1g
GLT-R 200	11.1	3.3	200
GLT-R 240	11.6	3.5	240
GLT-R 260	11.9	3.6	260
GLT-R 300	12.4	3.8	300
GLT-R 330	13.0	3.8	330
GLT-R 380	13.6	4.1	380



* Not available in Europe.



GLT-A Anatomical Gluteal Implant				
Product Code	Height (mm) ±2	Width (mm) ±2	Projection (mm) ±2	Volume (g) ±1g
GLT-A 215	13.6	99	2.8	215
GLT-A 285	13.9	10.9	3.3	287
GLT-A 305	14.3	11.0	3.4	305
GLT-A 345	14.5	11.4	3.6	345
GLT-A 385	15.0	11.7	4.2	385
GLT-A 440	15.2	12.1	4.6	440
GLT-A 470	15.2	12.1	4.8	470



* Not available in Europe.

Calf Implants

by EUROSILICONE®

- ✓ 2 configurations with 11 variations.
- ✓ Symmetric or asymmetric shape.
- ✓ High Cohesive Gel.
- ✓ Smooth shell surface.
- ✓ Disc on the back facilitates implant insertion.



ES 54		Symmetric Shape		
Product Code	Length (cm)	Width (cm)	Projection (cm)	Volume (cc)
54060	16.8	3.3	2.0	60
54090	17.0	3.7	2.8	90
54110	17.7	5.8	2.1	110
54120	18.9	5.2	2.6	120
54160	19.5	6.9	2.3	160
54180	19.0	5.9	3.1	180
54220	18.7	7.1	3.2	220




ES 54		Asymmetric Shape		
Product Code	Length (cm)	Width (cm)	Projection (cm)	Volume (cc)
54135	22.0	4.9	2.3	135
54140	18.6	5.3	2.6	140
54150	13.9	6.4	2.9	150
54165	23.9	5.6	2.3	165

Testicular Implants

by EUROSILICONE®

- ✓ 5 different variations.
- ✓ Oval shape.
- ✓ High Cohesive Gel.
- ✓ Smooth shell surface.
- ✓ Fixation tab eases suturing to surrounding tissues.




ES 520		Smooth Surface	
Product Code	Diameter (cm)	Length (cm)	Volume (cc)
52006	2.2	2.8	6.0
52012	2.8	3.3	12.0
52018	3.0	4.2	18.0
52022	3.3	4.7	22.0
52026	3.3	5.2	26.0

Testicular Implants

by **NAGÔR®**



- ✓ Filler options: Soft and semi – hard.
- ✓ Size options.
- ✓ Low profile reinforced suture loop.
- ✓ High performance silicones.
- ✓ 360° Barrier Shell Technology.



TFI		Firm		
Product Code	Size	Width (cm)	Length (cm)	Volume (g)
TFI-1	Small	2.6	3.1	15
TFI-2	Medium	3.4	3.9	27
TFI-3	Large	3.6	4.2	33


TGID		Gel-Fill		
Product Code	Size	Width (cm)	Length (cm)	Volume (g)
TGID-1	Small	2.6	3.1	13
TGID-2	Medium	3.4	3.9	23.5
TGID-3	Large	3.6	4.2	30

Facial Implants

Malar Implants

by **EUROSILICONE®**

- ✓ Medical-Quality silicone elastomer.
- ✓ Fixation holes to ease suturing to surrounding tissue.
- ✓ Smooth surface.



ES 33			
Product Code	Length (cm)	Width (cm)	Thickness (cm)
33001	2.4	1.7	0.3
33002	3.0	2.1	0.4
33003	3.4	2.3	0.4
33004	3.8	2.4	0.5

Nasal Implants

by EUROSILICONE®

- ✓ Multiple sizes available.
- ✓ Multiple perforations.
- ✓ High performance silicone elastomer.



ES 55			
Product Code	Length (cm)	Width (cm)	Base (cm)
55001	5.1	0.8	2.5
55002	5.3	0.8	2.8
55003	5.4	0.8	3.3
55004	5.9	0.9	3.5

Chin Implants

by EUROSILICONE®

- ✓ 2 configurations with 7 variations.
- ✓ Medical-Quality silicone elastomer.
- ✓ Smooth surface.



ES 32 / 34 / 36 / 38			
Product Code	Length (cm)	Width (cm)	Thickness (cm)
32000	3.4	1.3	0.9
34000	3.9	1.3	0.9
36000	4.2	1.5	0.9
38000	4.5	1.3	0.7



ES 325				
Product Code	Length (cm)	Depth (cm)	Width (cm)	Thickness (cm)
32501	6.2	3.0	0.6	0.6
32502	6.5	3.4	0.8	0.8
32503	7.0	3.8	0.9	0.8

Wound Management

Silgel™


62

Silgel™

Clearly better

The Silgel family of products has been developed to extend the concept of gel therapy by utilising high molecular weight silicone gel (polysiloxane) in the management of scars up to two years old.

Silgel and Silgel sheets are easy to apply and help to reduce and soften hypertrophic and keloid scars and associated erythema.



- ✓ Visibly reduces the dimensions of scars (length, width and height).¹
- ✓ Changes the consistency of the scar to produce softer tissue.¹
- ✓ Fades the red colouration of scars.¹
- ✓ Over 90% patients saw an improvement in the appearance of their scar within 16 weeks.¹

Silgel™ STC (Silicone Topical Gel)

Product Code	Box Quantity	Volume
STC-SE	20 Tubes	20ml Tube

Silgel™ Sheet Procedure Specific

Product Code	Box Quantity	Dimensions (cm)
TGS-AB (Abdominal)	1 Sheet	46.0 × 8.5
TGS-M (Sub Mammary)	2 Sheets	25.0 × 15.0
TGS-C (Circular)	2 Sheets	5.5

Silgel™ Sheet Square

Product Code	Box Quantity	Dimensions (cm)
TGS 1010	10 Sheets	10.0 × 10.0
TGS 2020	2 Sheets	20.0 × 20.0
TGS 4040	1 Sheet	40.0 × 40.0

Silgel™ Sheet Rectangular

Product Code	Box Quantity	Dimensions (cm)
TGS 1005	20 Sheets	10.0 × 5.0
TGS 1510	5 Sheets	15.0 × 10.0
TGS 3005	4 Sheets	30.0 × 5.0
TGS 3010	2 Sheets	30.0 × 10.0

¹Stewart, S.A., Dougall, G.M., Tafuro, E.M. The Use of Silgel STC-SE, a Topical Silicone Gel for the Treatment and Reduction of Hypertrophic and Keloid Scars. Plast Reconstr Surg Glob Open 2016; Vol 4: Issue 12: e1183

GCA Comfort Guarantee.



GCA Comfort Guarantee.



We offer a comprehensive and automatic warranty on all Nagor and Eurosilicone breast implants*

*On all silicone gel breast implants implanted by an appropriately qualified and licensed physician in accordance with accepted surgical procedures and techniques, and in strict accordance with GCA instructions for use as per product package insert.

What products are covered?

Silicone gel-filled breast implants.

What brands are covered?

Nagor and Eurosilicone.

Are previously implanted products covered?

Yes, the warranty covers all currently implanted Nagor and Eurosilicone gel-filled breast implants.

What does the warranty cover?

Provides replacement implants free of charge in the event of rupture of one or both implants.**

Provides replacement implants free of charge in the event of capsular contracture (Baker Grade III or IV) in one or both implants.**

** To learn more about the policy, please consult with your surgeon.

Replacements available

Any model, size or style of Nagor or Eurosilicone breast implant from the current catalogue.

Excess?

None.

How long does the guarantee last?

For the lifetime of the patient.

What information is required and why?

We need to verify that a Nagor or Eurosilicone product has been implanted and understand the nature of the claim. In order to do this we need information such as the model, size and lot number of the product implanted.

How long does it take to confirm if eligible for the warranty?

If the relevant information is provided the claim can be verified almost instantly.

Does the warranty cover rupture at any stage of surgery?

Yes.

Does the guarantee cover costs of surgery?

No, the guarantee is only on the product, it does not cover any surgical costs.

Do ruptured products need to be returned?

Yes, Competent Authorities require an inspection and investigation to be undertaken to establish the cause of the rupture.

Do intact products need to be returned?

No, please see the information below which highlights the steps that need to be followed for explanted intact products.

To Distributors, Hospitals, Clinics and Surgeons

Explanted intact product

It has been agreed that Nagor & Eurosilicone Mammary Implants which have been explanted, and which are intact and free of defects do not have to be physically returned to Nagor or Eurosilicone for examination.

How does it work?

1. A statement is provided from the Distributor/Hospital/Clinic/Surgeon that the device is intact and free of defects
2. Photographic evidence is provided of the device
3. A Customer Claim/Complaint Form is completed for each event

On receipt of the above, the complaint will be processed as normal and the guarantee confirmed, if applicable.

Dimensions & Volumes Chart.

Round Collection™
by EUROSILICONE®


IMPLEO™
by NAGÔR®




Round Collection™

by EUROSILICONE™

Volume g/cc for The Round Collection™




Depending on the product specification, 1g is approximately equivalent to 1cc.



by **NAGÔR®**

Abstract



For more information about our products
please visit: **www.gcaesthetics.com**
or email **info@gcaesthetics.com**

GLOBAL CONSOLIDATED AESTHETICS LIMITED.

Suite 601, Q_House, Furze Road, Sandyford Industrial Estate,
Dublin 18, Ireland.

EUROSILICONE®

EUROSILICONE® S.A.S.

Z.I de la Peyrolière
BP.68 84402 APT Cedex,
France.

T. +33 (0) 4 90 04 30 30

F. +33 (0) 4 90 04 60 06

NAGÔR®

NAGOR®

129 Deerdykes View, Westfield
Industrial Estate, Cumbernauld,
Glasgow. G68 9HN, UK.

T. +44 (0) 1236 780780

F. +44 (0) 1236 781234

f /gcaestheticsglobal

t /G_C_Aesthetics

in /company/gc-aesthetics

